www.bllifesciences.com

AngiX™
Interventional Systems

At the heart of critical moments

Interventional Systems



COMPONENTS

Manifold

Choice of 2, 3 or 4 Core
Available with color handle
Clear material for easy detection of bubbles
Large bore for maximum flow 500 PSI

Variants available on request.

Ordering Information

	On Handle	Off Handle
4 Core	B-MN-4RN	B-MN-4RF
3 Core	B-MN-3RN	B-MN-3RF
2 Core	B-MN-2RN	B-MN-2RF



Control Syringe

Available in 6ml, 10ml, 12ml & 20 ml
Low friction hub
Bold marking for correct measurement of dosage
Finger grips for one hand operation
Rotating male luer lock

Ordering Information

6ml	B-CSR-6RL
10ml	B-CSR-10RL
12ml	B-CSR-12RL
20 ml	B-CSR-20RL



Pressure Monitoring Line

PM Lines:
Choice of length from 6" to 84"
Choice of colour - Red, Blue, Yellow & Transparent
Minimum Priming Volume
ID 1.27 mm x OD 2.79 mm

Disposable Transducer Systems:
For Invasive Pressure monitoring
compatible with all Hemodynamic systems



Contrast Line

Available in 60" - 100"

Ordering Information

B-OT-26

High Pressure Line

Available in length 10", 20" or 48"
Braided PU material,1200 PSI rating
Rotating male luer lock
Also available in PVC material upon request

Ordering Information

10"	B-HPT10
20"	B-HPT20
48"	B-HPT48



Guide Wires

Variety of guide wires available in following variants:

O.D.	:	0.018", 0.021", 0.025", 0.032", 0.035", 0.038"
Length	:	45cm, 60cm, 70cm, 150cm
Material	:	Stainless steel wires with PTFE coated and Non coated, Nickle titanium
Tip	:	J tip, Straight



Ordering Information

For guide wire - Length 150cm, PTFE coated

	0.018"	0.021"	0.025"	0.032"	0.035"	0.038"
J Tip	B-GW018-150J	B-GW021-150J	B-GW025-150J	B-GW032-150J	B-GW035-150J	B-GW038-150J
Straight Tip	B-GW018-150ST	B-GW021-150ST	B-GW025-150ST	B-GW032-150ST	B-GW035-150ST	B-GW038-150ST

The Hydrophilic wires also available upon request. PTFE Guide wire available in 260 cm length upon request.

Introducer Needle / Puncture Needle

Highly precise, atraumatic bevel tip needle

Ordering Information

18G x 7cm	B-IN-18G7
20G x 4cm	B-IN-20G4
21G x 7cm	B-IN-21G7



Hemostasis Y:

Large bore
Compatible with catheters from 4Fr upto 9Fr and guide wires from 0.014" to 0.038" dia
Rated pressure 480 PSI
Spcl. feature of Click Y
Septum can be used in OPEN, SEMI OPEN and FULL CLOSED position

Options Available

Large Bore Rotating (Standard)
Y Click Type
Push / Pull

Customised Extension Lengths Available Upon Request

Ordering Information

Hemostasis Y Rotating (Standard)	-	BHYL - 03
Hemostasis Y Click Type	-	BHYL - 04
Hemostasis Y Push/Pull	-	BHYL - 05
Hemostasis Y Click Type includes Ext. Line 25cms - 4 way	-	BHYL - 04-10





Torque Device

Perfect grip for easy manipulation of guidewire

Compatible with Guidewire 0.009 to 0.0018 / 0.0018-0.0038

Ordering Information

BTD-02




Insertion Tool

Insertion Tool/Blunt Needle

Ordering Information

20GX10cm B-BN20G10



Inflation Device

Rated pressure 30 bar/430 PSI Ergonomic Handle

Available variants - Gun type / Straight

Ordering Information

B-ID-30



Radial Band

Ordering Information

BOT-92



CATHLAB KIT / ANGIO KIT

Pre-assembled kits for the convenience & quick set during cathlab procedures

Ordering Information/ Contents

BKT-53RN	Manifold 3 Port R/On, IV Set, Contrast Line, 10ml Control Syringe, 40" PM Line, 'Rotating' Haemostasic Y, Torque Device, Blunt Needle, 20" M/F extension & Injector Connector.
BKT-54RN	Manifold 2 Port R/On, IV Set, 10ml LL Syringe, 40" & 20" PM Line, 10ml syringes
BKT-342	Hemostasic Y 'Click-Press' type with 10" extension, Torque Device & Blunt Needle.
BKT-344	Hemostasic Y, 'Rotating ' type, Torque Device, & Blunt Needle
BKT-365	Accessories Kit Includes Haemostasic Y Click Type, Torque Device, Needle.
BKT-509	Haemostasic Y, 'Push-Pull' type with 10" extension, Torque Device, & Blunt Needle
BKT-506	Inflation Device 30atm, Haemostatic Y 'Rotating' type, Torque Device, & Blunt Needle
BKT-510	Accessories Kit includes Push/Pull Y with Extension 25 cms and 4 way Stopcock.



(Many Customised kits available upon request)

All the components are available individually packed or in assembled kit.

CATH LAB PROCEDURE PACKS

Customized Cath lab Procedure Packs are available as per user defined requirements and specifications.

Choice of wide range of components, customized pre-assembled cath lab kits and additional procedure accessories such as:

Towel clamp, Forceps, Swab gauze, Gallipot, Bowls (250ml, 450ml, 1000ml), Kidney Bowls (500ml, 700ml), Scalpel Blade, Sponge stick/swabs, Scissor, Fluro Cover, Trolley covers, Gowns (absorbent, non absorbent), Hand Towels, Cath Lab Drape with fenestration.



Customized Drape Set With Penetrations & Release Material
Contact Our Representative For Customized Drawings

INTRODUCER SETS FOR FEMORAL & RADIAL ACCESS

Size	Colour	Sheath			Dilator		Needle		Guidewire		Part No
		OD mm	ID mm	Length cm	Size guage	Length cm	Size guage	Length cm	Size inch	Length cm	
10 F	FUCHSIA	4.29	3.40	10	10F	15.0	18	7.0	0.035	45	BIS-10F18G35
9 F	BLACK	3.96	3.07	10	9F	15.0	18	7.0	0.035	45	B-IS-9F18G35
8 F	BLUE	3.58	2.74	10	8 F	15.0	18	7.0	0.035	45	B-IS-8F18G35
7 F	ORANGE	3.25	2.41	10	7 F	15.0	18	7.0	0.035	45	B-IS-7F18G35
6 F	GREEN	2.94	2.15	10	6 F	15.0	18	7.0	0.035	45	B-IS-6F18G35
5 F	GRAY	2.64	1.85	10	5 F	15.0	18	7.0	0.035	45	B-IS-5F18G35
5 F	GRAY	2.64	1.85	10	5 F	15.0	21	4.0	0.018	45	B-IS-5F21G18
4 F	VIOLET	2.28	1.49	10	4 F	15.0	21	4.0	0.018	45	B-IS-4F21G18
SIZES AND SPECIFICATIONS- RADIAL WITH HYDROPHILIC WIRE											
6 F	GREEN	2.83	2.20	7.4	6 F	13.0	20	4.0	0.025	45	B-IS-6F20GC25SSR*
5 F	GRAY	2.54	1.89	7.4	5 F	13.0	20	4.0	0.025	45	B-IS-5F20GC25SSR*

Seldinger technique.
Dilator locking prevents retrograde movement.
All sets includes syringe of 2.5 ml & surgical blade.

Variants available upon request.



B L Lifesciences Private Limited
A-245, Okhla Industrial Area Phase I,
New Delhi - 110020 INDIA.

+91 11 4510 0100
+91 11 4170 8464
blls@bllifesciences.com