
www.bllifesciences.com

VenX™
Central Venous Catheters

At the heart of critical moments


Central Venous Catheres Kits



VenX™
Central Venous Catheters

OUTSTANDING PERFORMANCE PROFILE FOR AN EFFICIENT TREATMENT

Specially Designed Soft & Conical Tip ensures Atraumatic, Safe & Easy Insertion





Biocompatible PU material is stiff & non kinking at room temperature, but softens at body temperature.




Radiopaque catheter & length markings ensures accurate placement.



Full range of Adult, Pediatric & Infant CVC Kits are available.



Components of Standard Kit



Kink resistant Nickel-Titanium Guidewire with soft & flexible J - tip offers better torque which helps in easy insertion & prevents vessel perforation.
Wires are Double Distal.
2 piece design of guidewire advancer.



Nickel-Titanium Guidewire

CVC kits are available with options of std. Needle, Needle with Back Flow valve & Introducer Syringe.



Puncture Needle Y Needle Introducer Syringe

VenX™
Central Venous Catheters

Ordering Information

	Catheter Size Ga/Fr	Lumen Ga	Cather Length cms	Dilator Fr	Guide Wire Nickle Titanium Double Distal	with Std. Needle	with Y Needle Std. Syringe	with Std. Needle Introducer Syringe
Single Lumen Catheter						P/N	P/N	P/N
1 Lumen 22Ga 4cms	22Ga	22	4	4	0.018x30 cm	B-CVC-22G01-04	---	---
1 Lumen 22Ga 4cms with 22G Cannula	22Ga	22	4	4	0.018x45 cm	B-CVC-22G01-04	---	---
1 Lumen 22Ga 8cms	22GA	22	8	4	0.018x45 cm	B-CVC-22G01-08	---	---
1 Lumen 22Ga 10cms	22Ga	22	10	4	0.018x45 cm	B-CVC-22G01-10	---	---
1 Lumen 22Ga 13cms	22 Ga	22	13	4	0.018x45 cm	B-CVC-22G01-13	---	---
1 Lumen 22Ga 40cms	22 GA	22	40	4	0.018x45 cm	B-CVC-22G01-40	---	---
1 Lumen 20Ga 13cms	20 Ga	20	13	4	0.018x45 cm	B-CVC-20G01-13	---	---
1 Lumen 18Ga 10cms	18 Ga	18	10	5	0.018x45 cm	B-CVC-18G01-10	B-CVC-18G01-10YN	B-CVC-18G01-10IS
1 Lumen 18Ga 13cms	18GA	18	13	5	0.018x45 cm	B-CVC-18G01-13	B-CVC-18G01-13YN	B-CVC-18G01-13IS
1 Lumen 18Ga 20cms	18 Ga	18	20	5	0.018x45 cm	B-CVC-18G01-20	B-CVC-18G01-20YN	B-CVC-18G01-20IS
1 Lumen 16Ga 13cms	16Ga	16	13	5	0.035x45 cm	B-CVC-16G01-13	B-CVC-16G01-13YN	B-CVC-16G01-13IS
1 Lumen 16Ga 16cms	16 Ga	16	16	5	0.035x45 cm	B-CVC-16G01-16	B-CVC-16G01-16YN	B-CVC-16G01-16IS
1 Lumen 16Ga 20cms	16 Ga	16	20	5	0.035x60 cm	B-CVC-16G01-20	B-CVC-16G01-20YN	B-CVC-16G01-20IS
1 Lumen 14Ga 20cms	14 Ga	14	20	8	0.035x60 cm	B-CVC-14G01-20	B-CVC-14G01-20YN	B-CVC-14G01-20IS
1 Lumen 16Ga 40cms	16GA	16	40	5	0.035x60 cm	B-CVC-16G01-40	B-CVC-16G01-40YN	B-CVC-16G01-40IS
1 Lumen 16Ga 60cms	16GA	16	60	5	0.035x70 cm	B-CVC-16G01-60	B-CVC-16G01-60YN	B-CVC-16G01-60IS
Double Lumen Catheter								
2 Lumen 4Fr 6cms	4 Fr	20, 22	6	5	0.018x 45 cm	B-CVC-402-06	B-CVC-402-06YN	B-CVC-402-06IS
2 Lumen 4Fr 8cms	4 Fr	20, 22	8	5	0.018x 45 cm	B-CVC-402-08	B-CVC-402-08YN	B-CVC-402-08IS
2 Lumen 4Fr 13cms	4 Fr	20, 22	13	5	0.018x 45 cm	B-CVC-402-13	B-CVC-402-13YN	B-CVC-402-13IS
2 Lumen 5Fr 8cms	5 Fr	18, 20	8	6	0.018x 45 cm	B-CVC-502-08	B-CVC-502-08YN	B-CVC-502-08IS
2 Lumen 5Fr 13cms	5 Fr	18, 20	13	6	0.018x 45 cm	B-CVC-502-13	B-CVC-502-13YN	B-CVC-502-13IS
2 Lumen 7Fr 13cms	7 Fr	14, 18	13	8	0.035x 45 cm	B-CVC-702-13	B-CVC-702-13YN	B-CVC-702-13IS
2 Lumen 7Fr 16cms	7 Fr	14, 18	16	8	0.035x 45 cm	B-CVC-702-16	B-CVC-702-16YN	B-CVC-702-16IS
2 Lumen 7Fr 20cms	7 Fr	14, 18	20	8	0.035x 60 cm	B-CVC-702-20	B-CVC-702-20YN	B-CVC-702-20IS
Triple Lumen Catheter								
3 Lumen 4.5Fr 6cms with 22G Cannula	4.5 Fr	20, 22, 22	6	5	0.018x45 cm	B-CVC-22C4.503-06	---	---
3 Lumen 4.5Fr 6cms	4.5 Fr	20, 22, 22	6	5	0.018x 45 cm	B-CVC-4.503-06	B-CVC-4.503-06YN	B-CVC-4.503-06IS
3 Lumen 4.5Fr 10cms	4.5 Fr	20, 22, 22	10	5	0.018x 45 cm	B-CVC-4.503-10	B-CVC-4.503-10YN	B-CVC-4.503-10IS
3 Lumen 5.5Fr 6cms	5.5 Fr	20, 22, 22	6	6	0.018x 45 cm	B-CVC-5.503-06	B-CVC-5.503-06YN	B-CVC-5.503-06IS
3 Lumen 5.5Fr 8cms	5.5 Fr	20, 22, 22	8	6	0.018x 45 cm	B-CVC-5.503-08	B-CVC-5.503-08YN	B-CVC-5.503-08IS
3 Lumen 5.5Fr 13cms	5.5 Fr	20, 22, 22	13	6	0.018x 45 cm	B-CVC-5.503-13	B-CVC-5.503-13YN	B-CVC-5.503-13IS
3 Lumen 7Fr 10cms	7 Fr	16, 18, 18	10	8	0.035x 45 cm	B-CVC-703-10	B-CVC-703-10YN	B-CVC-703-10IS
3 Lumen 7Fr 13cms	7 Fr	16, 18, 18	13	8	0.035x 45 cm	B-CVC-703-13	B-CVC-703-13YN	B-CVC-703-13IS
3 Lumen 7Fr 16cms	7 Fr	16, 18, 18	16	8	0.035x 45 cm	B-CVC-703-16	B-CVC-703-16YN	B-CVC-703-16IS
3 Lumen 7Fr 20cms	7 Fr	16, 18, 18	20	8	0.035x 60 cm	B-CVC-703-20	B-CVC-703-20YN	B-CVC-703-20IS
Four Lumen Catheter								
4 Lumen 8.5Fr 10cms	8.5 Fr	14,16,18,18	10	9	0.035x 45 cm	B-CVC-8.504-10	B-CVC-8.504-10YN	B-CVC-8.504-10IS
4 Lumen 8.5Fr 13cms	8.5 Fr	14,16,18,18	13	9	0.035x 45 cm	B-CVC-8.504-13	B-CVC-8.504-13YN	B-CVC-8.504-13IS
4 Lumen 8.5Fr 16cms	8.5 Fr	14,16,18,18	16	9	0.035x 45 cm	B-CVC-8.504-16	B-CVC-8.504-16YN	B-CVC-8.504-16IS
4 Lumen 8.5Fr 20cms	8.5 Fr	14,16,18,18	20	9	0.035x 60 cm	B-CVC-8.504-20	B-CVC-8.504-20YN	B-CVC-8.504-20IS

Option of Needle-less valve is also available in place of Injection port

B L Lifesciences Private Limited
A-245, Okhla Industrial Area Phase I,
New Delhi - 110020 INDIA.

+91 11 4510 0100
+91 11 4170 8464
blls@bllifesciences.com

CE Certified Products
ISO 13485

BL/CAT/04 Rev 02 Date 02-06-2025